
# SmartMat+ User Manual

## **Assistance Information**

For technical help, please call Podimetrics, Inc.

**Phone:** 1-888-498-6489

**Email:** [patients@podimetrics.com](mailto:patients@podimetrics.com)

**Mail:**

Podimetrics, Inc.

100 Dover Street

Somerville, MA 02144

Complaints related to device functionality, performance, or service can be submitted by phone or email.

## **About the Podimetrics SmartMat+**

The Podimetrics SmartMat+ includes sensitive electronics to measure weight, balance, and temperature over the soles of the feet as well as transmit that data wirelessly to a remote server for evaluation. Flexible temperature sensors are positioned directly under the cover and on top of a layer of foam that compresses when you stand on it. Each thermal scan that is taken with the SmartMat+ is transmitted wirelessly to a server where it can be reviewed by a healthcare provider for signs of inflammation on your feet.

The SmartMat+ can also provide information about weight and balance, measures that can be indicators of overall health.

# TABLE OF CONTENTS

<b>SAFETY INFORMATION</b> .....	<b>2</b>
<b>DESCRIPTION OF PARTS</b> .....	<b>4</b>
<hr/>	
<b>SETUP INSTRUCTIONS</b> .....	<b>6</b>
FIRST TIME SETUP SCREENS .....	8
PERSONALIZING YOUR RESULTS .....	10
<hr/>	
<b>DAILY USAGE INSTRUCTIONS</b> .....	<b>12</b>
DAILY SCANNING SCREENS .....	14
<hr/>	
<b>USING THE CONTROLLER</b> .....	<b>16</b>
MAIN MENU .....	17
DAILY SUMMARY .....	19
VIEW SCAN HISTORY .....	20
SCAN STATUS GUIDE .....	21
WATCH PERIOD .....	22
MEASURING BALANCE .....	24
OPTIONS .....	26
ABOUT DEVICE .....	27
<hr/>	
<b>DATA TRANSMISSION</b> .....	<b>28</b>
<b>CLEANING</b> .....	<b>29</b>
<b>BATTERIES &amp; CHARGING</b> .....	<b>30</b>
<b>WHAT TO DO WITH YOUR SMARTMAT+</b> .....	<b>31</b>
<b>INDICATOR LIGHTS</b> .....	<b>32</b>
<b>TROUBLESHOOTING</b> .....	<b>34</b>
<b>TECHNICAL INFORMATION</b> .....	<b>36</b>
<hr/>	
<b>MANUAL DE USUARIO</b> .....	<b>39</b>

# **SAFETY INFORMATION**

## **Indications For Use**

The Podimetrics SmartMat+ is part of the Podimetrics Remote Temperature Monitoring System™ (RTM System). The RTM System is intended to be used by a patient in conjunction with a healthcare professional or caretaker for periodic evaluation of the temperature over the soles of the feet for signs of inflammation.

The RTM system provides information indicating when the patient and healthcare provider should communicate for further evaluation and treatment regarding any persistent localized inflammation observed on the feet via the electronic sensing system and remote visualization of its data.

The Podimetrics RTM System is intended to be used in the home under the direction of a healthcare professional as an adjunct to, and not in replacement of, self-examination and periodic foot care and examination conducted by a healthcare professional and does not diagnose any specific disease state. It is for prescription use only.

The SmartMat+ also measures weight and balance, which can be indicators of overall health. Weight and balance data is for informational purposes only.

## **Contraindications**

Patients with an active ulcer or other sores that have broken the skin on the soles of their foot or feet should not scan the wounded foot on the Podimetrics SmartMat+. Patients may continue using the Podimetrics SmartMat+ to monitor an unwounded foot.

## **Essential Performance**

The essential performance of the SmartMat+ is to initiate a scan when pressure is applied to the FootMat; display scan progress; reliably record foot temperature, weight, and balance data; and transmit the data to a remote server.

## Benefits

It can be difficult for patients at risk for inflammatory foot diseases to evaluate the soles of their feet for signs of inflammation associated with an injury. The temperature data transmitted wirelessly by the Podimetrics SmartMat+ may be used by you and your doctor to identify signs of inflammation in your feet. With this information, your doctor may choose to initiate a physical examination and actions to avoid complications like foot ulcers, infection, gangrene, and amputation.

## Risks

The Podimetrics SmartMat+ does not cause any pain or discomfort to use. The top layer has been tested for skin biocompatibility. If damaged or used improperly, inaccurate temperature readings may be taken, which could affect the clinician's evaluation of inflammation.

## Warnings

- Do not use the SmartMat+ if it appears damaged or contaminated.
- Do not use the SmartMat+ for more than one person because it may lead to cross-contamination or user data mix-up.
- If any irritation, redness, or rash develops on the bottom of your feet, discontinue use and contact Podimetrics immediately.
- No modification of this equipment is allowed.

## Precautions

Do not place other objects or coverings on the SmartMat+ when it is being used. Do not wear socks, slippers, or shoes while scanning on the SmartMat+. The mat's measurements may not be accurate if used with an external covering or while wearing socks or shoes - scans must be taken while barefoot.

# DESCRIPTION OF PARTS

---


## SmartMat+ (PMHS300U)

The SmartMat+ includes a **FootMat**, **Ramp** (detachable), **Controller** and a **Charger** (for the Controller).

## FootMat (PMHS400U)

The FootMat measures your foot temperature, weight, and balance when stepped on for scanning.

The Ramp smooths transitions from floor to FootMat and is detachable.


**⚠ CAUTION:** Do not wear socks or shoes when scanning your feet on SmartMat+.

---

## Controller (PMHS300U)

The Controller communicates wirelessly with the FootMat to display scans and scan results. It can be used while attached or detached from the FootMat.

 **Note:** The Controller is not a touchscreen. Use the buttons to navigate.

### Navigation Buttons

**OK** - Make and confirm selections.

**Back** - Exit or go back to previous screens.

**Up / Down** - Move up or down through menu options, lists, or information screens.

### Indicator Lights

See pg 32 for Indicator Lights.

### Charger

The Charger is for recharging the Controller.

The FootMat uses AA batteries and can't be recharged.

(See pg 30 for more on charging or batteries.)


# SETUP INSTRUCTIONS

## 1 Remove clear protective covers

Remove any clear protective covers from the FootMat and Controller. They will interfere with scanning if not removed.

## 2 Power ON the Controller

Use the switch on the back to turn the Controller ON.

**Note:** Do not plug the Controller in at this point. It comes pre-charged.

## 3 Power ON the FootMat

Look at the bottom of the FootMat


- (1) If there is a plastic tab sticking out of the battery compartment, pull it all the way out and throw it away.
- (2) Then use the switch to turn the FootMat ON.

## 4 Position the FootMat

Set the FootMat on a flat floor where you won't trip on it and where you can comfortably stand for 20 seconds.

The dark blue side should be facing up, and the ramp should be facing you.


**CAUTION:** Do not wear socks or shoes when scanning your feet on SmartMat+.


**Note:** Device is ready for use just seconds after switching ON.

## 5 Position the Controller

The Controller can be attached at the top of the FootMat, or it can be detached and placed on a nearby table or shelf.

*Attached to FootMat*

*Example: Detached from FootMat & placed on a shelf*

## 6 Follow the steps on the Controller screen

The Controller screen will walk you through the rest of the setup process (see pg 8, First Time Setup for more detail).

## 7 Leave devices powered ON

Leave the Controller and FootMat powered ON at all times so they can send your scans and receive scan results.

### **⚠ CAUTION:**


- Do not place SmartMat+ in a walkway like a hall or doorway where you might trip on it.
- Do not place SmartMat+ in a location where it is likely to get wet, like next to a shower.
- Do not place SmartMat+ in direct sunlight or next to heat sources like vents or space heaters.
- Place SmartMat+ in a location that minimizes the risk of tampering by children or pets.

**🔗 Note:** Keep the box for storing the SmartMat+.

# FIRST TIME SETUP SCREENS


---

The following screens are displayed during the setup process. Make sure the FootMat is turned on during the setup process.


## 1 Welcome

Software will load when the Controller is powered on.

## 2 Connecting to Cell Network



SmartMat+ will automatically try to connect to a cellular network.

## 3 Scan Your Feet


Take your first foot scan.

- Place your bare foot or feet on the FootMat to start the scan.
- Stay still while the scan is in progress (~20 seconds).
- “Scan Saved” will appear when the scan is done.

## 4 Weekly Scan Log

Shows which days of the week you've completed a scan.


## 5 Weight

Shows your weight after each scan.


**⚠ CAUTION:** Only stand on SmartMat+ without using assistance if you feel steady and comfortable enough doing so.


**6 Weight Change**  
Shows your change in weight since your last scan.


**7 Balance**  
Shows if you're Steady or Unsteady compared to your personalized Stability Baseline (see pg 24).

**8 Battery**  
Shows the Controller's battery life.


**9 Personalize Your Scan Display Now?**  
Choose if you want to personalize your scan display now.

**Personalize Your Scan Display**  
Choose which measures you'd like to see or hide after each scan.

- Use the UP or DOWN buttons to move between lines.
  - Press OK to select Show or Hide.
  - Press OK or "Done" to continue.
- You can also change these later in the Options menu (see pg 26).

**10 Goodbye**  
The screen will turn itself off. Leave the Controller and FootMat powered ON at all times so they can upload scans and receive scan results.

**Note:** Accurate weight and balance measurements require standing on the SmartMat+ for the entire scan (~20 seconds) without using any support. If you are unable to stand on the SmartMat+ without assistance, you may wish to hide the Weight and Balance displays during setup or in the Options menu (see pg 26).

# PERSONALIZING YOUR RESULTS


---

## Foot Template


Before you can see your daily foot scans, Podimetrics needs to create a personalized foot template that matches your feet shape and size.

Creating your foot template may take up to 8 scans, so please scan daily.

While your foot template is being created, all scans will show as Results Pending. After your template is completed, scans will be marked as Results Pending until your scans are processed (typically around the end of the next business day), then their statuses will change to indicate the results. Please continue to scan daily to see your progress.

*Results Pending in Daily Summary*


*Results Pending in Foot Scans*

---

## Steadiness Baseline

In order to see information about your balance, Podometrics will create your personalized Steadiness Baseline, which may take up to 8 scans. The SmartMat+ will not display information about your balance until your Steadiness Baseline has been created. Once your Steadiness Baseline is complete, the SmartMat+ can let you know if you are more or less steady than your baseline after each scan (see pg 24).

**Note:** You can choose to show or hide your balance results at any time by changing Balance Display in the Options menu (see pg 26). The SmartMat+ always collects balance data so that it is available to view at any time.

# DAILY USAGE INSTRUCTIONS

---

**Note:** For best results and most reliable weight and balance measurements, follow this process at the same time every day.

## 1 Perform daily foot exam

Examine your feet to ensure there are no signs of developing ulcers or open wounds.

**WARNING:** If you notice broken skin or a developing ulcer, stop using the SmartMat+ and contact Podometrics immediately.

## 2 Check the FootMat

Before use, inspect the FootMat for wetness, damage, or discoloration. Make sure that the FootMat is dry before using it. If it is damaged, do not use the SmartMat+ and contact Podometrics for a replacement. If it is soiled or discolored, try to clean off any residue using a damp cloth. If the discoloration cannot be removed, contact Podometrics to see if a replacement is necessary.

## 3 Stand barefoot on the FootMat

Scan your bare foot/feet on the SmartMat+. Scanning starts automatically when you step on the FootMat.

For best results, scan with your eyes open and your bare feet facing forward in a relaxed position and placed in the middle of the blue sections of fabric. If you're scanning a single foot, you can place your foot on either the right or left side of the blue fabric. Your feet should always be facing forward (toes pointing away from the ramp).

**CAUTION:** Do not use with water or lotion on your feet.

## 4 Complete the scan

Stay as relaxed as possible and keep your feet still for 20 seconds while the scan completes. The Controller will show a progress bar and notify you when the scan is done.

**Note:** Although use while standing is preferred, the SmartMat+ can also be used while sitting, including while in a wheelchair, if enough pressure is applied to complete a scan. Do not use the SmartMat+ in a way that you feel may increase your risk of falling or tripping.

**Note:** Keep your foot/feet still during the entire scan. Moving your feet during the scan can interrupt the scan. If this occurs, step on the FootMat again to restart a scan.

**Note:** While taking a scan, you can hold the Controller, place it on a table or shelf, or attach it to the FootMat.

**Note:** Weight and balance measurements will not be accurate when taken while sitting to scan or while holding onto anything for assistance. If you sit to scan or always hold onto something like a wall, table, or walker while standing on the SmartMat+, you may wish to hide the Weight and Balance displays in Options (see pg 26).

## 5 Leave Controller and FootMat powered ON after use

After completing your daily scan, leave both the Controller and FootMat powered ON. There is no need to turn the devices off – they will go to sleep after each scan.

**Note:** Keep the Controller near the FootMat at all times to ensure a strong connection and so that you can watch the screen while scanning.

# DAILY SCANNING SCREENS

To start a scan, stand barefoot on the FootMat and a scan will begin with the screens below.

**⚠ CAUTION:** Do not wear socks or shoes when scanning your feet on SmartMat+.


You don't need to hold or use the Controller to start scanning. It will automatically wake up when you step on the FootMat.


**1 Welcome**

**2 Stay Still**  
A reminder to stay still for the scan.


**3 Scan Progress**  
Shows the progress of the 20 second scan.

**4 Weekly Scan Log**  
Shows which days of the week you've completed a scan.

**Note:** Units displayed on the device and in this guide reflect the user's settings (lb or kg). Examples shown may use pounds.


- 5 Weight (optional)**  
Shows your weight.


- 6 Weight Change (optional)**  
Shows your change in weight since your last scan.

- 7 Balance (optional)**  
Shows if you are Steady or Unsteady compared to your personalized Stability Baseline (see pg 24).

- 8 Battery**  
Shows the Controller's battery life.

- 9 Goodbye**  
The screen will turn itself off. Leave the Controller and FootMat powered ON at all times so they can upload scans and receive scan results.

 **Note:** Optional displays can be set to Show or Hide during First Time Setup or in Options (see pg 26).

# USING THE CONTROLLER


---

## Waking Up the Controller

Pressing any button on the Controller will it wake up and display key information, such as reminders, important messages or your daily scan results.


You can use the buttons to go to the Main Menu or set the Controller down and it will go back to sleep in a short while.

### Screens that may appear when the Controller wakes up:


#### Reminders

Displays when you step on the FootMat.


#### Important Messages

Alerts you when you're in a Watch Period.

#### Scan Results

Displays your daily scan results.

#### Main Menu

Access to other functions, such as Daily Summary, View Scan History, Options or About Device.

 **Note:** When you're done using the Controller, do not power it off. You can set it down and the screen will automatically go to sleep after a short period of inactivity.

Leave the Controller and FootMat powered ON at all times so they can upload scans and receive scan results.


# MAIN MENU

---

The Main Menu screen provides key information at a glance.


Battery  
Level

Cellular  
Signal

## Recent Scans

Shows the number of scans you've taken that day and if you have any results pending.

## Reminders or Watch Period Alert

If you have a Reminder or are in a Watch Period, they will appear here instead of the Recent Scans information (see pg 22 for Watch Period).

The Main Menu also provides access to additional areas and functions. Use the buttons to navigate the Main Menu and other screens:


- **UP / DOWN** - to move up or down through menu options, lists, or information screens.
- **OK** - to make and confirm selections.
- **Back** - to exit or go back to previous screens.

 **Note:** Ensure the SmartMat+ has a cellular signal so that it can send your scans and receive scan results.

# MAIN MENU (continued)


---

You can access these additional areas and functions from the Main Menu. Each area is further described on the next pages.


## Daily Summary

See a summary of each day's scan result, Weight, Weight Change and Balance.


## View Scan History

View your daily scan results, weight trends, balance trends, and a calendar view highlighting days you scanned.

## Options

Change settings for: Daily Summary displays (i.e. Weight, Weight Change, Balance), reminder light timing, display sleep timer, language or units.

## About Device


Access to device information, manual software update, factory reset, or view device logs.

# DAILY SUMMARY


---

The Daily Summary screen provides a snapshot of the metrics for each day. You can press UP or DOWN to navigate the dates.

The optional fields can be set to Show or Hide during First Time Setup or in Options (see pg 26).

*Optional fields set to Show*


*Optional fields set to Hide*

- 1 Scan Status:** Shows if you took a scan that day.
- 2 Scan Result:** Shows the status of that day's scan (see pg 21 for Scan Status Guide).
- 3 Weight (optional):** Shows your weight for that day.
- 4 Weight Change (optional):** Shows your change in weight between that day's scan and the previous scan.
- 5 Balance (optional):** Shows if you were Steady or Unsteady compared to your personalized Stability Baseline for that day (see pg 24).

 **Note:** If optional fields are set to Hide, SmartMat+ will still log them in case you want to view them later.

# VIEW SCAN HISTORY


In View Scan History you can review your foot scan results, weight trend, balance and scan log. You can press UP or DOWN to navigate the dates for each.

## Foot Scans


View a history of your Foot Scans.

Selecting an entry and pressing OK will bring up an image of that day's scan and a scan result explanation.


## Weight Trend

View your weight trends over 7, 14, 21, or 28 days, as well as yearly.


## Balance

View your balance compared to your baseline over time.

## Scan Log

View the days you successfully scanned. These are days marked by a checkmark or are highlighted in green.


# SCAN STATUS GUIDE

---


Scan Image	Icon	Description
	 Results Pending	Scan received but not processed yet
	 Normal Scan	Scan doesn't show areas of high temperature
	 Elevated Scan	Scan shows temperatures that are different than normal
	 Unclear Scan	Scan shape didn't match your personal template
No Image	 Scan Taken	Scan was taken and is being uploaded for processing
No Image	 No Scan Taken	No scan taken
No Image	 No Data	No data available

# WATCH PERIOD

---

*Red Watch Period alert*

*Red Watch Period banner*

Watch Period is triggered if you have two consecutive Elevated Scans (dual feet scanners) or a single Elevated Scan (single foot scanners).

Because your scans show a pattern of high temperature, you may be at higher risk of developing a foot wound for the next 21 days. Pay extra attention to your feet during a Watch Period, continue to scan daily, and talk to Podimetrics to learn how to prevent sores or wounds on your feet. If you continue to have Elevated Scans, your Watch Period may be extended beyond 21 days.

The SmartMat+ measures the temperatures on the bottom of your feet/foot. Research shows that higher foot temperatures may mean a sore is starting. It could also be a sign of an infection or another problem.

We recommend you do the following:

- Talk to Podimetrics Patient Support to learn how to prevent a sore or wound on your foot.
- Continue to scan daily. We will continue to watch your scans every day to see if the temperatures stay high. Check your feet every day for swelling, calluses, blisters, color changes, or wounds. If you notice any of these, notify Podimetrics Patient Support Immediately. Ask someone for help with checking your feet if you can't see well.

- 
- Always wear your prescribed shoes, both inside and outside.
  - Decrease activity.
  - Elevate your feet when sitting down.

By checking the temperature of your feet regularly, we can find these changes and help prevent bigger problems.

 **Note:** The system cannot detect or diagnose foot ulcers, open wounds, or any other diseases.

# MEASURING BALANCE

---

The SmartMat+ contains highly sensitive force transducers that measure small shifts in your body position as you stand on it and complete your scan. By measuring these shifts in body position, the SmartMat+ calculates your average Center of Pressure velocity (COPv), or how fast your body sways as it tries to stay balanced while standing still. Higher COP velocities suggest greater instability and an increased risk of falling.

Every time you scan, SmartMat+ can tell you whether you're more or less steady than your "normal" stability. To start, SmartMat+ uses your first 4-8 scans to create your personalized Stability Baseline - think of this as defining your "normal" stability. Once your baseline is complete, SmartMat+ can let you know if you are more or less steady compared to your personalized Stability Baseline every time you scan.

If SmartMat+ shows a message like "Unsteady Compared to Your Baseline" or "Take Extra Care Today", that's just a helpful nudge to go a little slower and be more mindful as you move around that day. It's not something to panic about, it's a tool to help you stay aware and confident as you move through your day.

 **Note:** Accurate stability measurements require standing on the SmartMat+ for the entire scan (~20 seconds) without using any support. If you are unable to stand on the SmartMat+ without assistance, you may wish to turn off the Balance displays during device setup or in the Options menu (see pg 26).

 **CAUTION:** Only stand on the SmartMat+ without using assistance if you feel steady and comfortable enough doing so.

---

To ensure the most reliable stability measurements, it is recommended you do the following:

- Scan at the same time every day.
- Scan the same way every day.
- Stand with your bare feet comfortably positioned on the FootMat.
- Use a quiet, relaxed stance with eyes open, arms at your sides and a slight bend in your knees.
- Do not try to keep a rigid pose; just stay relaxed.

### Balance

We're building your baseline  
This may take 8 scans  
Continue scanning daily

### Building your baseline

Your baseline is being created.  
This will take 4-8 scans.

### Steady

compared to your baseline

### Steady

Your measured stability matches your baseline.

### Unsteady

compared to your baseline

### Unsteady

Your measured stability is less than your baseline. **Stay safe during your day:**

- Use your cane, walker, or any support you normally rely on.
- Wear safe, non-slip shoes.
- Take things a little slower. Avoid rushing or carrying too much at once.

Take extra care today

 **Note:** Balance data is for informational purposes only.

# OPTIONS

---

In Options, you can customize the settings below.

For all settings: Press OK to toggle or confirm a selection.

## Daily Summary


Choose which optional displays you want to Show or Hide following your scans and in your Daily Summary views.

 **Note:** The device will still collect data for optional displays that are set to Hide. If you switch them back to Show, you'll be able to see the logged history.


## Daily Reminder Light

Set your preferences for the Daily Reminder Light: On or Off; adjust its brightness; or when it should turn on (Wake Time) or off (Sleep Time).


## Standby Timer

Set how long the screen stays on after a period of inactivity.

## Language

Change the default language shown on the Controller. (Not all languages may be available.)

## Units

Choose to display results in pounds or kilograms.

# ABOUT DEVICE


---

In About Device you can view details about the system to help with troubleshooting or with support calls.

## Device Info


View the current software and build version of your SmartMat+.

## Check for Updates


Manually check for any updates.

 **Note:** SmartMat+ will automatically check for and install updates regularly.

## Reset to Factory Settings

Resets SmartMat+ to factory settings, erasing all user data and preferences. To do this, you will need to call Podmetrics to obtain a PIN.

## Device Logs

Shows device activity logs, which can help during support calls.

# DATA TRANSMISSION

---

 **Note:** WiFi is not needed. Data is uploaded via cellular networks securely and automatically.

## Data uploads automatically

If you have placed SmartMat+ in a location with a good cellular signal, SmartMat+ will automatically upload data to Podometrics. Cellular coverage may not be available in all regions. Contact Podometrics if you have difficulty finding a location with a good signal. Data will not be uploaded outside of the United States. The SmartMat+ will save at least 31 scans in memory if outside of cell service, and it will retry uploading scans daily.

# CLEANING

---

## FootMat Cleaning

Clean the FootMat as needed. Wipe the top surface of the FootMat with lukewarm water and a clean cloth to remove any dirt or residue that may build up and affect the sensor reading or make it slippery. Let the FootMat completely air dry before use. If it is difficult for you to clean the FootMat while it is on the floor, you can place it on a table.

## Controller Cleaning

Clean the Controller as needed. Wipe the surface and buttons of the Controller with lukewarm water and a clean cloth to remove any dirt or residue.

**⚠ CAUTION:** Do not use chemical cleaners on the FootMat or Controller as it may damage them.

**⚠ CAUTION:** Do not clean the FootMat with soap as it may leave a slippery residue. Do not immerse the FootMat or Controller in water.


**⚠ CAUTION:** If liquid spills on the FootMat or Controller, wipe and clean surfaces immediately with a clean cloth.

# BATTERIES & CHARGING

---

## Charging the Controller

The battery level of the Controller is shown after each use. To charge it, plug the included charger into a standard wall outlet and insert the USB-C end of the cable into the charger port on the side of the Controller. It takes up to four hours to fully charge the Controller.

### **⚠ CAUTION:**

- Do not connect any device or cord other than those provided by Podimetrics to the USB-C port. Using other cables may cause damage.
- Make sure the charging cord is not in a place where someone could trip on it.
- Keep the charging cable out of reach of children. The charging cable may be a strangulation risk.

## FootMat Batteries

The FootMat does not require charging – it comes with batteries pre-installed. If the Controller alerts you that the FootMat batteries are empty, you can replace the four AA batteries by following the instructions below:

1. Turn the FootMat over and locate the battery door.
2. Use a small phillips head screwdriver to loosen the screw holding the battery door.
3. Replace all four AA batteries with fresh batteries – ensure the batteries are correctly oriented (+ and -).
4. Close the battery door and tighten the screw back.

# WHAT TO DO WITH YOUR SMARTMAT+

---

## Traveling

Feel free to take your SmartMat+ with you to use it at your destination. You can detach the Ramp from the FootMat to make it easier to pack.

**Especially when flying, turn the FootMat and Controller OFF while in transit.** If you don't plan to bring your SmartMat+ with you, call Podimetrics to inform them that you won't be using it. The SmartMat+ is not designed for use outside the United States.

**⚠ CAUTION:** The charger may not work in other countries.

## Storage

Store your SmartMat+ powered off in its shipping box in a cool, dry place out of direct sunlight. The materials may degrade over time if exposed to sunlight or high temperatures. After cold (-4 °F/-20 °C) or hot (113 °F/45 °C) storage, allow SmartMat+ to reach room temperature (68 °F/20 °C) for 15 min before use.

**⚠ CAUTION:** When not using the FootMat, keep it flat on the ground and where it won't be a trip hazard; do not store it leaned against a wall or in an upright position to prevent it from falling.


## Power Isolation

To isolate the device (fully disconnect from power), unplug the charger from the wall outlet.

**⚠ CAUTION:** Only use a wall outlet you can easily reach to unplug when needed.

# INDICATOR LIGHTS


If the information here does not resolve the issue or if you experience any malfunction or injury, please call Podometrics Support at 1-888-498-6489.

## Controller Light Guide

Light	Meaning	Action
<b>A</b>  <b>Blinking Green</b>	Controller battery is charging.	No action required.
<b>A</b>  <b>Solid Green</b>	Controller battery is fully charged.	No action required.
<b>A</b>  <b>Solid Yellow</b>	Controller is plugged in but not charging because of an error or overheating.	<ol style="list-style-type: none"> <li>1. Unplug the Controller.</li> <li>2. Wait 10 minutes.</li> <li>3. Plug the Controller back in.</li> <li>4. If problem persists, contact support.</li> </ol>
<b>B</b>  <b>Blinking Blue</b>	Reminder to scan.	Take your daily scan.
<b>B</b>  <b>Solid Purple</b>	Controller is initializing.	No action required.

## Indicator Light

## FootMat Light Guide

Indicator Light	Meaning	Action
 <b>Blinking Blue</b>	FootMat is looking for the Controller.	No action required.
 <b>Solid Yellow</b>	FootMat is not connected to the Controller, OR Footmat is not working properly.	<ol style="list-style-type: none"> <li>1. Place FootMat and Controller close together.</li> <li>2. Turn both OFF, wait 10 seconds, then turn both ON.</li> <li>3. If problem persists, contact support.</li> </ol>
 <b>Blinking Green</b>	Your scan is in progress.	No action required.
 <b>Solid Green</b>	Device is ready and working.	No action required.
 <b>Solid White</b>	Scan completed, going back to sleep.	No action required.
 <b>No light</b>	FootMat is asleep, OR Footmat is not working properly.	<ol style="list-style-type: none"> <li>1. Step on FootMat to wake it up.</li> <li>2. If there is still no light, contact support.</li> </ol>

# TROUBLESHOOTING

---

If the information here does not resolve the issue or if you experience any malfunction or injury, please call Podometrics Support at 1-888-498-6489.

Code	Problem	Meaning	Action
none	Screen shows “Stay Still”	Movement was detected preventing the scan from starting or completing.	Stay still for the scan to start and also stay still for the entire scan (~20 seconds).
none	Screen shows “Apply More Pressure”	The scan cannot start or complete because there’s not enough pressure detected.	Apply more pressure when stepping on the FootMat.
1	The Controller display won’t load	Possible issue with firmware.	Switch the Controller OFF and ON to restart it.  Check if Controller is charged.
2	Network connection fail	A connection couldn’t be established.	Try moving the system to a new spot, such as near a window.
3	Unable to connect FootMat	The Controller cannot communicate with the FootMat.	Switch the FootMat OFF and ON, or the FootMat may need new batteries.
4	Scan could not start	Possible insufficient pressure detected.	If you’re sitting for the scan, try applying more pressure.

<b>Code</b>	<b>Problem</b>	<b>Meaning</b>	<b>Action</b>
<b>5</b>	<b>Scan interrupted</b>	Movement while scanning will interrupt the scan.	Try to stay still for the entire scanning. Step on the FootMat to restart a scan.
<b>6</b>	<b>Scan save failure</b>	Saving was interrupted.	Step on the FootMat to restart a scan.
<b>7</b>	<b>Cellular SIM not active</b>	Issue with cellular SIM.	Contact support.
<b>8</b>	<b>Support required</b>	An issue is present that requires contacting support.	Contact support.
<b>9</b>	<b>SmartMat+ deactivated, unable to scan</b>	Your SmartMat+ is deactivated.	Contact support.
<b>10</b>	<b>Weak cell signal</b>	Your SmartMat+ has not been connected to a cellular signal for an extended amount of time. The system can't continue to perform without a cell signal.	Try moving the system to a new spot, such as near a window.

# TECHNICAL INFORMATION

---

## Safety Symbols

Do not use the product if the packaging is damaged.

Read the operating instructions carefully.

**IP21**

FootMat is rated IP21. It is protected against inserted fingers and dripping water.

**IP22**

Controller is rated IP22. It is protected against inserted fingers and dripping water.

The SmartMat+ is a Type BF applied part and provides a high degree of protection against electric shock.

This product contains electronics.  
Do not dispose of in municipal waste facilities.

Rx only

This product is for prescription use only.

Device uses RF transmitter that produces non-ionizing radiation.

Product is single patient reusable.

Storage and transport -20°C to 45°C.

---

## Specifications

The SmartMat+ is designed for use in the home in ordinary indoor conditions. It must not be used outdoors, in wet areas, or in proximity to flammable anesthetic mixtures with air or oxygen/nitrous oxide.

**Operating Range:** 59 to 104 °F (15 to 40 °C)

**Storage and Transport:** -4 to 113 °F (-20 to 45 °C)

**Operating, Storage and Transportation Humidity:** 15 to 93% RH, noncondensing

**Operating Pressure:** 700 hPa to 1060 hPa atmospheric pressure

**The Maximum Enclosure (Accessible Surface) Temperature:** 106 °F (41.7 °C)

---

**Temperature Accuracy:** ± 1 °F (± 0.6 °C)

**Resolution:** 0.2 °F (0.1 °C)

**Normal Scan Time:** 20 seconds

**Spatial Resolution of Sensors:** 0.56 inches (~14 mm) vertical and horizontal spacing

---

**Maximum User Weight:** 400 lbs (182 kg)

**Weight Accuracy:** ± 1.0 lbs (± 0.45 Kg)

**Weight Resolution:** 0.1 lbs (0.05 kg)

---

**Device Dimensions:** FootMat w/ramp: 18.5 x 15 x 1 inches  
Controller: 6.3 x 3.9 x 1.4 inches

**Device Weight:** FootMat w/ramp: 8 lbs (3.6 kg)  
Controller: 0.6 lbs (0.3 kg)

**Service Life:** 3 years of normal, daily use

---

**LCD Display:** 4.3 inches LCD

---

**Cellular Radio:** LTE Cat-M1 radio, 700-1980 MHz,  
200mW maximum output power  
FCCID: XMR2021BG950AGL

**BlueTooth Radio:** Bluetooth Low Energy (BLE) radio, 2.4 GHz,  
1.66mW maximum output power  
FCCID: Y82-DA14531MOD

# TECHNICAL INFORMATION (continued)

---

**Controller Battery:** Rechargeable lithium ion battery pack,  
3.4 to 4.2 V, 1850 mAh,  
About 14 days of daily use

**FootMat Battery:** 4 replaceable AA Alkaline Batteries  
2 yrs of daily use

**Controller Battery** Medical Power Supply  
**Charger:** Input Rating: 100-240VAC, 50-60 Hz  
Output Rating: 10W, 5.0V@2.0A  
Mfg. P/N: GlobTek WR9HA2000USBCFMEDR6B  
Cable length: 60 inches

---

**Product Certifications:** IEC 60601-1 (2005) +A2 (2020), IEC 60601-1-2:2014,  
IEC 60601-1-6 (2010) +A2 (2020),  
IEC 60601-1-11 (2015) +A1 (2020), FCC Part 15  
Subpart B Class B and Canada ICES-003 Class B

---

**Patents:** <https://www.podimetrics.com/intellectual-property>

## Guidance and Manufacturer's Declarations and Recommended Separation Distances

See LBL-0096 Rev A - Supplemental technical information sheet.

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